



GI UK  Endo Therapy

POLYP TRAPS

OUR RANGE OF POLYP TRAPS ARE DESIGNED TO SIMPLIFY RETRIEVAL & IDENTIFICATION OF POLYPS

P-TRAP

Polyp traps are integral to collection of small lesion pathology. Losing or damaging tissue samples has a detrimental impact on patients & procedures. GI UK Medical have a range of polyp traps to suit all clinical requirements.

SOFT CATCH TRAYS

The soft catch tray system permits effective retrieval of pathology without compromising suction and loss of procedure time. The perforated drawers reduce risk of biomaterial exposure. The snap-lock system ensures drawers are securely closed to the trap, reducing risk of drawer handling error.



MAGNIFYING GLASS DESIGN

Large draw system & magnifying glass design, allows clear identification of polyp retrieval.



- Quick and easy set up to aid theatre process flow.
- Integrity of pathology maintained with soft catch drawer design

Product Code	Description	Minimum Delivery Quantity
VTL/PT-02	2 Drawer Polyp Trap	20

5-CHAMBER POLYP TRAP

PRODUCT INFO

- 5-chamber rotatable trap permits multiple samples to be collected without loss of suction.
- 5 clearly marked chambers to match polyp excision location.
- Transparent chambers enables identification of polyp.
- Intuitive design allows simple and fast set up to support theatre flow



Product Code	Description	Minimum Delivery Quantity
BP-40702	5 Chamber Polyp Trap	10
PT-01-05	5 Chamber Polyp Trap	10
VTL/PT-05	5 Chamber Polyp Trap	25

SINGLE CHAMBER POLYP TRAP

- Single chamber trap offers an economic alternative for cases where limited pathology is to be collected.
- The clear transparent bubble case permits unobstructed view to confirm sample has been collected.



Product Code	Description	Minimum Delivery Quantity
BP-40070008	1 Chamber Polyp Trap	10
PT-01-01	1 Chamber Polyp Trap	10
VTL/PT-01	1 Chamber Polyp Trap	50



 Endo Therapy

 Endo Decontamination

 Centurion PPE

 Capsule Endoscopy

 Watchman Managed Inventory

 Endo Essentials

GET IN TOUCH

 43 Park Place, Leeds. LS1 2RY

 0113 513 4870

 info@giukmedical.com