

Retrieval Forceps



Description	Grasper Jaw Length	Grasper Max. Opening	Catheter Ø	Catheter Length	Minimum Ø of Working Channel	Box Quantity	Product Code
Rat Tooth Jaw	5.1mm	6mm	2.3mm	180cm	2.8mm	5	CJ-PPDR-23-180
Alligator Jaw	5.8mm	9mm	2.3mm	230cm	2.8mm	5	CJ-2418GH

Rat-Tooth Alligator Grasping Forceps

Outer Catheter Ø	Working Length	Long Cups	Box Quantity	Product Code
FOR BRONCHOSCOPES WITH WORKING CHANNELS FROM 2.0MM				
1.8mm	120cm	5.8mm	10	BP-40090021
1.8mm	120cm	5.8mm	10	GFR-ALI-18-120-SU
FOR GASTROSCOPES AND COLONOSCOPES WITH WORKING CHANNELS FROM 2.8MM				
2.3mm	160cm	7.12mm	10	BP-40090022
2.3mm	230cm	7.12mm	10	BP-40090023
2.3mm	160cm	7.12mm	10	GFR-ALI-23-160-SU
2.3mm	230cm	7.12mm	10	GFR-ALI-23-230-SU

