

Type II Medical Face Masks



Type IIR Medical Grade Surgical Face Mask (Non-Sterile)

- Disposable Type IIR EN14683 Face Mask
- High Bacterial Filtration Efficiency (BFE) 98%
- Manufactured to EN14683 type II standards
- Protective 3 Layer Construction
- Ear Loop, Pleat Style Mask with integrated nose clip
- Pack of 50 Single Use Masks

Colour	Size	Type	Adjustable Nose Clips	Box Quantity	Product code
Blue	Universal (17cm x 9cm)	3-layer Ear Loop Mask	Yes	50	TYPEII*

